

X-DRS VET Expert

Technical Specifications



Generator

Power	30 kW
Anode	Rotating, max. current 425 mA
Focus	0,6/1,3 mm
Power supply	230 V; 16 A (3 kW)
Ma-range in steps of 25	25 to 300
KV range in steps of 1KV	40 to 120
Collimator	R104

X-ray system

Floating table-top	
Dimensions with table top (W x H x D)	63 x 25 x 31 inch
Dimensions with stand (W x H x D)	90 x 79 x 53 inch
Weight of the entire system	250 kg

Detector

Detector	X-DR XL
Pixel pitch	139 µm
Pixel matrix	3,072px x 3,072px
Weight	4.0 kg
Safety class	IPX1
Surface load	150 kg
Point load	100 kg to Ø 40 mm

Software

X-AQS PACS VET

X-AQS PACS Vet
X-AQS Module Archiving on CD/DVD
X-AQS Module Data Compression
includes for client
X-AQS Module Patient CD/DVD
X-AQS Module Import Bitmap Images
X-AQS Module Import DICOM Images
X-AQS Module Import AQS Images
X-AQS Module Image Export
X-AQS Module Mail Export

X-AQS Client

1 workstation (expandable up to 4 workstations)
X-AQS Viewer
X-AQS Diagnosis
License per Client, floating license
for X-AQS Software