

# X-R Portable 20 Plus

## Technical Specifications



High-Frequency X-Ray Generator	
Max. output	1.6 kW @ 80 kV
Output	100 kV / 20 mA
kV range	40 - 100 kV (1 kV steps)
mAs range	0.3 - 50 mAs (22 steps)
Current output range	12 - 20 mA
Tube type	Toshiba D-124
Focal spot	1.2 mm x 1.2 mm
HU capacity	20 kHU
Cooling rate	250 HU/s
Collimator lamp	50 W Halogen
Collimator auto timer	30 s
SID scale length	2 m
Laser pointer	Dual laser
Power input	220 / 240 V
Power frequency	50 / 60 Hz, single phase
PROM memory	4 APR
Dimensions (L x W x H)	352 mm x 220 mm x 200 mm (13.9" x 8.7" x 7.9")
Weight	9.8 kg (21.6 lbs)
Included accessories	Aluminum transport case